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## Instruction for use of SELF-BREAKING SCREW SYSTEM

Self-breaking SWING OFF/ BRING OFF SCREW SYSTEM medical devices are made according to the harmonized standards identified in the essential safety 93/42/CEE (D.L. 46/97) e 2007/47/CE (D.L.37/2010).

**Caution:** only highly qualified orthopaedic surgeons are allowed to implant the SELF-BREAKING SCREW systems in suitably equipped places. The manufacturer and the distributors shall not be responsible for missed compliance to the specifications designated and for an incorrect use of the product.

### Description

SELF-BREAKING screws are designed to provide fixation of internal fractures of small dimension or to bone reconstruction to rapidly reach a correct osteointegration of the same. However, these devices are not intended to be load-bearing or to withstand excessive activity, therefore the patient shall limit the physical activities that could lead to excessive mechanical loads or to the displacement of the fractured edges delaying the healing. SELF-BREAKING screws are SINGLE-USE devices

**N.B.** For the dimensions, check the product label.

### Material

All materials used for implant products comply with the regulations for medical use:

Screws: Titanium Alloy Ti6Al4V ELI (ISO 5832-3)

### Indications

SELF-BREAKING screws are intended for fixation of fractures and bone reconstruction. Some examples are:

- Osteosynthesis of small bones
- Weil osteotomies
- monocortical fixation
- Osteotomies and fixation of foot and hand fractures

### Contraindications

SELF-BREAKING screws are contraindicated in the following cases:

- sensitivity to foreign bodies, hypersensitivity (proven or presumed) to the materials of the devices
- local or systemic inflammation acute or chronic
- active or latent infection
- severe osteoporosis
- systemic or metabolic pathologies
- physiological or anatomical anomalies

### Warnings and precautions

Surgeons who intend to use the SELF-BREAKING SCREW medical devices, should be aware of the following:

- Implants which are not load bearing structures are unable to support bone loading relating to normal physical activities. These devices are not intended to be the substitution of the normal bone structure. Therefore patients must be provided with adequate external support systems or limit physical activity which could cause excessive mechanical strain or the displacement of the device. For this reason post-surgery instructions are extremely important.
- The correct choice of the implant relating to type and size is essential for the success of the treatment.
- Knowledge of pre-operative and intra-operative procedures, the good fracture reduction and the correct screw placement are essential factors for the success of the implant. Surgeons are advised to follow all the phases indicated in the surgical technique.
- The components belonging to SELF-BREAKING SCREW system produced by Overmed shall be used in combination with the dedicated instruments according to the surgical technique prescribed by the Manufacturer.
- When choosing the patient the following factors which can jeopardize the correct post-operative progress must be considered: obesity, work or activities which entail stressful situations, alcoholism, drug abuse, smoking or degenerative diseases; reluctance to accept the possibility to perform multiple surgical intervention for revision or substitution; infectious diseases and compromised healing. The surgeon is responsible of the selection criteria for the patients.
- The SELF-BREAKING SCREW medical devices are single-use. After removal or after implant attempts all components of the devices must be disposed of according to hospital procedures. Do not reuse already inserted implants. Overmed assumes no responsibility for such use. Reuse of disposable devices may cause infection or threaten the device functional performance.
- The surgeon should be discuss with the patient, before the surgical intervention, all possible risks, precautions, advices, consequences, complications and adverse reaction due to the surgical procedure and device implant. The surgeon shall give all the necessary information to the patient, before the intervention.

### Surgical Technique

A surgical protocol is available for the use of this device. It's the surgeon's responsibility to use the instruments recommended in the surgical technique prescribed by Overmed. Never attempt a surgical session with defective, damaged or suspect instruments. Inspect all components before use.

### Adverse Effects

Possible adverse effects are pain, discomfort, nerve or soft tissue damage due to the presence of an implant or due to surgical trauma. Fracture of the implant may occur due to excessive activity, prolonged loading upon the device, incomplete healing, or excessive force exerted on the implant during insertion. Implant migration and/or loosening may occur. Metal sensitivity, histological, allergic or adverse foreign body reaction resulting from implantation of a foreign material may occur.

### Sterilization and packaging

SELF-BREAKING SCREWS implants package includes:

- The decontaminated implant.
- N° 5 labels, check integrity.
- This instruction sheet.

### IMPORTANT

The SELF-BREAKING SCREWS are provided in STERILE packages (double coated medical packets). STERILE is indicated clearly on the label on the package. HEADLESS screws are sterilized (25 kGy beta irradiation) and sterilization has a five-year warranty.

Products indicating STERILE can be taken into the operating theatre and opened by specialized personnel, who must carefully open the external bag and reverse the internal packet containing the STERILE device without touching the sterile parts as indicated by the person in charge.

**Important:** A color change indicator label indicates the radiation exposure: check the indicator to assure that ray exposure has been carried out and verify the integrity of the double packaging before utilizing the product. In case of damage do not use the product and return it to Overmed. The manufacturer declines any responsibility for re-sterilized device by the user.

**Do NOT use SELF-BREAKING SCREWS already inserted.** Reuse of device may cause high risk of infection, contamination and implant failure.

### Implant identification

Identification of the implant is guaranteed by the information indicated on the label on the packaging. Apply the sticky label on the patient's case history folder, in order to assure traceability of the information in case of negative patient feedback.

### Damaged SELF-BREAKING SCREW medical devices


In the case of damage which could jeopardize the patients' health, once removed return the product to the supplier together with all identification labels (type of implant and code, manufacturer batch code, decontamination code) in order to activate immediately the complaint procedure and corrective action.

### Implant removal

The surgeon who implants the device is responsible for the evaluation of the clinical follow-up. As a general rule, all risks must be carefully considered and advantages and disadvantages relating to the removal or possibility to leave the device in the body, discussed with the patient.

**WARNING:** The self-breaking screws are not suitable for fixation with screws of posterior elements (spinal peduncle) of cervical, dorsal or lumbar spine

### Symbol legend

	Medical Device
	Product code number
	Screw diameter in millimeters
	Screw lengths in millimeters
	Material
	Batch code
	Use-by date (YEAR-MONTH)
	Do not re-use
	Caution
	Sterilized using beta rays
	Do not use if package is damaged
	Do not re-sterilize
	Manufacturer