FRA — Dispositif Médical: MY-SCREW
REF : 38.88.XX

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ENG — Medical Device: MY-SCREW
REF: 38.88.XX

13 sizes available - XX length from 10 mm to 34 mm
Material: titanium ISO 5832-3

Indications for use

MY-SCREW are designed to hold bone fragments in place during the bone healing process

MY-SCREW has been specifically designed for foot surgery, in particular for use in fixation of the chevron osteotomy for hallux valgus correction. The screw can be used in both cortical and cancellous bone.

Contra indications, Warning and possible adverse effects

Do not use in case of allergy to one of the material components.

General surgical contra indication: poor quality bone due to osteoporosis...

The implants should never be reused. It's essential that the surgeon is fully conversant with the prescribed surgical technique and that due consideration is given to the information provided below.

Improper patient selection, with particular regard to patient weight and functional demands, as well as placement and positioning of the implant may result in unusual stress conditions and subsequent reduction in their functional life.

The correct handling of implant is paramount. **Contouring or altering the device should be avoided.** Implants must not be modified or adapted.

Premature loading or muscular activity must be avoided as this may be problematic.

The surgeon should select the appropriate size and dimensions of the implant to suit the physiological functions of the patient and the type of osteosynthesis to be treated.

The product should NOT be implanted in the human body in contact with other metallic materials or different chemical composition as this may cause an electrolyte reaction and corrosion may occur.

The patient should be instructed to report any unusual changes around the operated site to his physician (e.g.: sign of infection) In case of non respect of those indications and contra indications, possible adverse effects included infection, subclinical nerve damage, arterial puncture, tendon risk, metal sensitivity, loosening, bending, or fracture of the device and/or the bone.

The patient must state that he/she is fitted with a metallic implant during MRI examinations or scans. The professional in charge of the examination shall ascertain the safety of the implant and patient, along with the risk of artefacts and data loss.

SURGICAL TECHNICAL

The specific ancillary equipment, especially the screwdriver blade (ref: 14.30.21) is required for fitting ORTHO CAPE.

1. Realization of the chevron osteotomy.
2. Establishment of the proximal pin. It should be bicortical in the diaphysis of the metatarsal. Oriented from proximal in distal to medial in lateral.
3. Move the head using a lever inserted into the first metatarsal
4. The proximal pin is then pushed into the head.
5. The second screw is parallel to the first
6. Exostosectomy.
7. Akin Osteotomy


Do not use if the packaging is damaged.

Device disposal: In the event of explantation, the implant should be given to a specialist department to ensure that it is disposed according to strict environmental and health regulations. A product explanted for deficiency should be decontaminated and returned to the manufacturer.

Sterilization method: gamma irradiation

 ORTHO CAPE 13 rue Mespoul,
31400 TOULOUSE - FRANCE.

 Products should not be stored for a long time at temperatures below 10 °C or above 33 °C.

 The devices are not intended to be sterilized or reused. Resterilization or re-use can cause infections or product performance loss. Preserve from humidity to protect the readability of the label.

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