## REUSABLE INSTRUMENTS FOR FOOT SURGERY

Reference: IFU-RIFS Version C – November 05, 2024

### English

#### Product Description/Intended Use

This instructions for use provides information on the care, cleaning, disinfection, maintenance and sterilization of reusable surgical instruments manufactured by Ortho Cape.

Reusable Instruments for Foot Surgery are provided non-sterile and must be cleaned, disinfected and sterilized prior to the first, and every subsequent use. For proper cleaning, disinfection and sterilization, follow the reprocessing instructions provided in this IFU.

The following cleaning and sterilization instructions should be integrated with the health care facility procedures. National rules, regulations and / or restrictions must be included in this process.

Instrument boxes do not provide a sterile barrier: sterilization packaging should be used to keep the device sterile.

The instruments are delivered non-sterile in transport containers. The transfer sheet (which specifies the non-sterile status of the surgical instrument) is supplied with each ancillary for incoming quality control and prior to sterilization.

Reusable Instruments should exclusively be used for the placement or explantation of ORTHO CAPE implants.

For specific instructions on surgical use, please refer to the Instructions for Use that accompany the corresponding Ortho Cape Implant, or contact your local Ortho Cape representative to obtain a surgical technique. Reprocessing instructions of Reusable Instruments for Foot Surgery manufactured by Ortho Cape are provided below.

#### Warnings

- The instructions provided in this guide have been validated by Ortho Cape. It is the responsibility of the healthcare facility to ensure that processing is performed using the required equipment, materials and personnel at a defined processing area. Equipment and processes must be validated and monitored regularly.
- To minimize the risk of infection, these devices must be cleaned and sterilized prior to the first use and after every subsequent use.
- Assembled instruments will not clean or sterilize completely. Disassemble all instruments prior to cleaning and sterilization.
- Using unspecified sterilization cycles or agents may damage the device or result in incomplete sterilization.
- Wear appropriate protective equipment individual as gloves, eye protection, etc.
- Reusable Instruments are intended to be used at ambient temperature. It is recommended to not use this material directly after steam sterilization and to wait for the return to ambient temperature.
- The following precautions should be observed by hospital staff working with contaminated or potentially contaminated medical devices. Pointed or sharp devices should be handled with utmost caution.
- Do not use metal brushes or scouring pads during manual cleaning procedures. Use soft nylon bristle brushes and cleaning brushes of suitable sizes for the devices to be processed.
- Highly alkaline cleaning detergents can create stains on anodized surfaces and loss of elasticity in the case of silicone parts.
- Cleaning detergents containing acid neutralizing agents or natron should not be used for cleaning sterilization containers.
- The use of demineralized water prevents the formation of stains and corrosion; for this purpose, we recommend the use of disinfectants with corrosion protection.
- Do not use cleaning agents containing iodine, iodide, active chlorine, chloride, mercury aldehyde or bromine, bromide: these products are corrosive.

- Do not use mineral oil or silicone lubricants as they coat microorganisms, prevent direct contact of the surface with the steam and are difficult to remove.
- Drills, pins, rasps and sharp instruments should be carefully inspected after treatment to ensure that the sharp edges are not damaged.
- Screwdriver blades should be carefully inspected after treatment to ensure they are not bent or dull.
- Some instruments are surgically invasive devices intended for transient use such as drills and taps. These instruments should be used with cautions. Rotation may give rise to conflict with the soft tissues or areas in contact.
- Pointed or sharp devices should be handled with utmost caution.
- The use of some instruments to be fitted on motors (drill, burr, ...) may give rise to potential overheating between the instrument and the bone. It is recommended to spray these instruments with physiological saline solution during use.
- Do not cross-sterilize the device. Using multiple sterilization methods may significantly reduce the performance of the device.
- Do not leave the device in solutions longer than necessary. This may accelerate normal product aging.
- Improper processing will negate the responsibility of Ortho Cape for the performance of the device.

#### Reprocessing

These reprocessing instructions are provided in accordance with AAMI TIR 12, AAMI TIR 30, ISO 17664-1 and ISO 17665-1. While they have been validated by the manufacturer of the medical device as being capable of preparing the device for re-use, it remains the responsibility of the processor to ensure that the reprocessing as actually performed, using equipment, materials, and personnel in the reprocessing facility, achieves the desired result. This normally requires validation and routine monitoring of the process. Ortho Cape recommends users observe these standards when reprocessing medical devices.

#### Pre-treatment

- Cleaning the devices within 30min after use and should be transported in closed or covered containers to avoid contamination;
- It is recommended for disassemblable instruments, (example: disassemble AO blade from screwdriver handle), to disassemble these devices for decontamination in order to achieve the desired level of sterility assurance.
- Ortho Cape surgical instruments can be cleaned with alkaline or enzymatic detergents, specific for materials for hospital use, considering the raw material of construction.

#### Handing Washing

1. Prepare alkaline or enzymatic detergent according to the manufacturer's recommendations.
2. Completely immerse the instruments in the enzymatic solution and leave to soak for 20 minutes +/- 1 min. Use a soft nylon bristle brush to gently clean the tool to eliminate all visible dirt and areas difficult to access, example: Articulated part, cavities, asperities.
3. Cannulated instruments must be freed of any bone residues that may be inside the central canal of the device. It is recommended to pass a k-wire inside the cannulated instruments before going through the ultrasonic bath.
4. Rinse with tap water for 3 min +/- 10 seconds all the device.

#### Ultrasonic Washing

1. After the operations referred to in handing washing, immerse the instruments in the ultrasonic bath.

## REUSABLE INSTRUMENTS FOR FOOT SURGERY

Reference: IFU-RIFS Version C – November 05, 2024

2. Leave to soak in the solution subjecting it for 10 min at 45-50 kHz in the ultrasonic bath. Use only suitable detergent with broad antimicrobial spectrum, neutral pH and compatible with raw material of the reusable instruments (stainless steel, titanium, and polypropylene).
3. Rinse instruments in purified water for at least 3 minutes or until all blood and bone remains have been removed.
4. Repeat the steps of cleaning the ultrasonic bath and rinse above and then remove the moisture with a clean, lint-free cloth

### Automated Washing

1. After the operations referred in the handwashing, remove the devices to rinse them in purified water for at least 1 min and then rinse the entire device
2. Wash the devices for 10 minutes at 93 ° C with a detergent solution suitable for washer / disinfectant
3. Rinse with deionized water and dry.

### Inspection before sterilization:

After cleaning it is necessary to visually check the presence or absence of contaminants; if so, the cleaning cycle must be repeated. The functional test and / or visual check must be performed after cleaning and assembling the previously disassembled parts. The instruments must be examined in order to highlight mechanical damage (e.g. fractures, deformations, corrosion, etc.) Damaged instruments must be discarded and replaced. If it presents corrosion, discoloration, staining or damage contact Ortho Cape to replace such defective device.

### Sterilization:

Sterilization is not a substitute for cleaning. Instruments must be thoroughly cleaned prior to sterilization ;

Sterilization parameters may change when using a sterilization tray. Refer to the Instructions for Use provided with the tray for appropriate parameters.

- Cleaned and dried instruments must undergo a validated sterilization cycle before use in the surgical field.
- Gather the device if it has been disassembled for decontamination.
- Steam sterilization packages or bags may be used.
- Instrument sets must be properly prepared and packed in trays and / or boxes to allow steam to penetrate and enter into direct contact with all surfaces.

### Steam sterilization

Reusable Instruments for Foot Surgery manufactured by Ortho Cape shall be sterilized by moist heat (steam). Time and temperature parameters required for sterilization vary according to the type of sterilizer, cycle design, and packaging material. Please review the sterilizer manufacturer's sterilization instructions, or hospital procedures, prior to sterilization. The minimum recommended parameters for steam sterilization are as follows:

	Temperature	Pressure	Exposure Time	Wrapping
Europe	134 °C	2.2 bar	18 min	Double Wrapped
USA	132°C (269.6°F)	2.2 bar	4 min	Double Wrapped
UK	134°C	2.2 bar	3 min	Double Wrapped

Note: A double protective packaging method type AAMI / ISO 11607-1/2 is recommended.

### Storage

When handling packaged cases, take care to avoid damaging the sterile barrier. The healthcare establishment should define a storage shelf-life for packaged instruments, on the basis of the sterile packaging and the recommendations of the manufacturer of this packaging.

### Hospital responsibilities with regards to Ortho Cape instruments loan / Number of potential reuse

Orthopedic surgical instruments generally have a long service life and may therefore be used multiple times. However, this service life may rapidly be reduced due to poor handling or unsuitable protection. Instruments that are no longer performing must be returned to Ortho Cape.

The instruments must undergo all stages of decontamination, cleaning, disinfection, inspection and final sterilization before being returned to Ortho Cape. The decontamination documentation must be provided with instruments returned to Ortho Cape.

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