FRA — Dispositif Médical : VIS EAGLE 7 MM / VIS ARTHROSTAB 7 MM

REF : VIS EAGLE 7.0 VX38.99.XX

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ENG — Medical Device : EAGLE 7 MM SCREW / ARTHROSTAB 7 MM SCREW

REF : EAGLE 7.0 SCREW VX38.99.XX

(A : thread 16 mm; B : thread 22 mm; C : thread 28 mm)

REF : ARTHROSTAB 7.0 SCREW VX38.22.XX

(XX length from 30 to 90 mm)

13 sizes available: XX length from 10 mm to 34 mm

Indications for use:

The use of these screws is mainly in 3 indications: Osteotomy of the calcaneus (+++ translation, subtraction, addition) / Substrate arthrodesis / Talo crural arthrodesis

Contra indications, Warning and possible adverse effects

Do not use in case of allergy to one of the material components

General surgical contra indication: poor quality bone due to osteoporosis. . .

The implants should never be reused. It's essential that the surgeon is fully conversant with the prescribed surgical technique and that due consideration is given to the information provided below.

Improper patient selection, with particular regard to patient weight and functional demands, as well as placement and positioning of the implant may result in unusual stress conditions and subsequent reduction in their functional life.

The correct handling of implant is paramount. Contouring or altering the device should be avoided. Implants must not be modified or adapted.

Premature loading or muscular activity must be avoided as this may be problematic.

The surgeon should select the appropriate size and dimensions of the implant to suit the physiological functions of the patient and the type of osteosynthesis to be treated.

The product should NOT be implanted in the human body in contact with other metallic materials or different chemical composition as this may cause an electrolyte reaction and corrosion may occur.

The patient should be instructed to report any unusual changes around the operated site to his physician (e.g.: sign of infection)

In case of non respect of those indications and contra indications, possible adverse effects included infection, subclinical nerve damage, arterial puncture, tendon risk, metal sensitivity, loosening, bending, or fracture of the device and/or the bone.

The patient must state that he/she is fitted with a metallic implant during MRI examinations or scans. The professional in charge of the examination shall ascertain the safety of the implant and patient, along with the risk of artefacts and data loss.

SURGICAL TECHNICAL

(The specific ancillary equipment is required for fitting ORTHO CAPE implants)

1. Osteotomy of the calcaneus

Patient in lateral decubitus, ankle to operate on top + pneumatic tourniquet.

Horizontal posterior incision. Dissection and protection of the sural nerve.

Transverse osteotomy of calcaneus under scopic control. Translation of the calcaneus (as medial, lateral or vertical case down or up).

EAGLE 7 screw / ARTHROSTAB 7 mm screw, kit spigot placed, under scopic control back and forth between the tuberosity and the calcaneus body, to maintain the reduction.

The measurement provided provides the measurement of the intraosseous path of the spindle.

We then perform the bone drilling about 2/3 of the length of the pin with the wick, taking care to cross the osteotomy focus under scopic control.

The stepped wick is then introduced to the stop to prepare the area of the screw head.

Lastly, the EAGLE 7 screw / ARTHROSTAB 7 mm screw is inserted corresponding to the spindle measurement minus 15 mm.

This EAGLE 7 screw / ARTHROSTAB 7 mm screw is introduced by controlling its progression in scopie, so that the 2 threads of the screw is well on both sides of osteotomy to stabilize and ensures its compression.

The spindle can then be removed.

Verification of the stability of the assembly as well as the compression with the scopic control.

Additional immobilization by a removable boot at 90 °.

2. Substrate arthrodesis

Patient in lateral decubitus, ankle to operate on top + pneumatic tourniquet

A lateral horizontal incision is made in front of the tip of the fibrilar malloleus allowing exposure of the subtalar joint.

After sectioning of the interosseous ligament, dis-assembly of the subtalar and removal of the cartilaginous articular surfaces with the chisel to be struck and the posterior and posterior talc Gouge forces are performed.

Preparation of the congruence of the surfaces and confrontation of the latter, under perioperative scopic control.

Placement of two pins (provided in the kit EAGLE 7 / ARTHROSTAB 7 mm) introduced in plantar under the heel by a short cutaneous approach.

These pins are placed, under scopic control, parallel to each other in anterior posterior, and direction up and forward, ascending. They will allow during screwing to compress the 2 surfaces prepared in a harmonious way.

The intraosseous path of the pins is then measured using the measuring bar provided in the EAGLE 7.0 mm / the ARTHROSTAB 7 mm kit (measurement by subtraction).

Drilling is then performed on approximately 2/3 of the length of the first posterior pin in bone area with the wick, taking care to cross the joint to be treated under scopic control.

The stepped wick is then introduced to the stop to prepare the area of the screw head.

Finally, the posterior EAGLE 7.0 mm screw / ARTHROSTAB 7 mm screw corresponding to the spindle measurement minus 15 mm is introduced.

This screw EAGLE 7.0 mm / ARTHROSTAB 7 mm screw is introduced by controlling its progression in scopie, so that the 2 threads of the screw is well on both sides of the joint to merge and ensures the compression of its surfaces.

3. Talo crural arthrodesis

Supine patient with pneumatic tourniquet

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Joint opening with protection of the vasculo-nervous bundle, visualization of the osteoarthritic tibial joint

Performing extensive arthrolisis to visualize joint surfaces. It is made with the chisel to strike and the Cauchoix a cartilaginous decrease of the articular surfaces.

The two medial and lateral aligners with release of the osteophytes are then released by reciprocating sawing.

A fibula shortening osteotomy corresponding to the cumulative resection area of the articular surfaces of the tibia and talus is then performed.

These osteotomies make it possible to face the articular parts perfectly.

Maintaining the arthrodesis in place and in good position by two previous compressive staples.

An anteroposterior talian catch-up screw introduced by the anterior cortex of the tibia is added, obliquely downwards and behind.

To do this, we introduce a spindle of the screw kit EAGLE 7 / ARTHROSTAB 7 mm, under scopic control along the axis defined before taking care not to enter the subtalar joint.

The measurement provided provides the measurement of the intraosseous path of the spindle.

We then perform the bone drilling about 2/3 of the length of the pin with the wick, taking care to cross the joint well to be treated under scopic control.

The stepped wick is then introduced to the stop to prepare the area of the screw head.

Lastly, the EAGLE 7 mm screw / ARTHROSTAB 7 mm screw is

inserted corresponding to the spindle measurement minus 15 mm.

This EAGLE 7 mm screw / ARTHROSTAB 7 mm screw is introduced by controlling its progression in scopie, so that the 2 threads of the screw is well apart and other of the joint to merge and ensures the compression of its surfaces.

The spindle can then be removed.


Verification of the stability of the assembly as well as the compression with the scopic control.

Registry in optional residual line spacing. Additional immobilization by a removable boot at 90 °

Do not use if the packaging is damaged

Device disposal: In the event of explantation, the implant should be given to a specialist department to ensure that it is disposed according to strict environmental and health regulations. A product explanted for deficiency should be decontaminated and returned to the manufacturer.

Sterilization method: gamma irradiation

Products should not be stored for a long time at temperatures below 10°C or above 33°C.

The devices are not intended to be sterilized or reused. Resterilization or re-use can cause infections or product performance loss. Preserve from humidity to protect the readability of the label

inserted corresponding to the spindle measurement minus 15 mm.

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