


BRING-OFF

Ø 2.0 ET 2.7 MM

- Vis en titane ISO 5832-3
- Sécable, auto perforante et auto taraudeuse
- Avant-pied
- Conditionnement stérile



Ø 2.0 mm	
LONG.	RÉF.
10	A1502010
11	A1502011
12	A1502012
13	A1502013
14	A1502014
15	A1502015

Ø 2.7 mm	
LONG.	RÉF.
10	A1502710
11	A1502711
12	A1502712
13	A1502713
14	A1502714
15	A1502715
16	A1502716
17	A1502717

ANCILLAIRE

N°	PRODUIT	RÉF.
1	Lame de tournevis 3 pans pour vis Bring-off	K1010035
2	Manche de Tournevis AO	14.32.03

