

MY SCREW

Ø 2.5 MM



Ø 2.5 mm

- Vis en titane ISO 5832-3
- Auto taraudeuse, auto foreuse et auto compressive
- Avant-pied
- Conditionnement stérile
- Broche Ø 1.0 pour vis Ø 2.5 mm

Ø 2.5 mm	
LONG.	RÉF.
12	38.88.12
14	38.88.14
16	38.88.16
18	38.88.18
20	38.88.20
22	38.88.22
24	38.88.24
26	38.88.26
28	38.88.28
30	38.88.30
32	38.88.32
34	38.88.34

ANCILLAIRE

N°	PRODUIT	RÉF.
1	Forêt canulé pour vis My Screw Ø2.5	14.37.25
2	Forêt canulé pour vis My Screw Ø2.5 XL	14.37.27
3	Manche de tournevis AO universel	14.32.03
4	Tournevis (lame) Torx 7	14.30.22
5	Mesureur sur broche 70 mm	14.33.22
6	Tournevis (lame) Torx 7 canulé	14.30.21
7	Broches Ø 1.0 Lg 70 mm	14.63.01
8	Boite de rangement	14.10.88

