

# SBS

Ø 2.0



Ø 2.0mm

- Vis en titane ISO 5832-3
- Sécable, auto perforante et auto taraudeuse
- Avant-pied
- Conditionnement stérile

Ø 2.0 mm

LONG.	RÉF.
09	14.52.09
11	14.52.11
12	14.52.12
13	14.52.13
14	14.52.14
15	14.52.15
17	14.52.17
19	14.52.19
21	14.52.21
24	14.52.24

## ANCILLAIRE

N°	PRODUIT	RÉF.
1	Lame de tournevis SBS	14.30.02
2	Manche de Tournevis AO	14.32.03

