




# Ortho CAPE

## FOOT SURGERY


For any questions regarding our products,  
contact us at [info@ortho-cape.com](mailto:info@ortho-cape.com)

# PERCUTANEOUS BURRS

OrthoCAPE




- X90CrMoV18 Steel burrs
- Percutaneous surgery
- Singly sterile packaging
- Standards burrs Ø 2.32 Lg 70



	N°	REF.	Ø	CUTTING LG	TOTAL LG
SHANNON SHORT	①	38.37.01	2 MM	8 MM	70 MM
SHANNON LONG	②	38.37.02	2 MM	12 MM	70 MM
SHANNON ISHAM	③	38.37.03	2 MM	12 MM	70 MM
SHANNON XL	④	38.37.07	2 MM	22 MM	70 MM
SHANNON XXL	⑤	38.37.08	3 MM	20 MM	70 MM
SHANNON XXL LG 30	⑥	38.37.09	3 MM	30 MM	100 MM
BROPHY	⑦	38.37.06	5 MM	15 MM	70 MM
WEDGE Ø3,1	⑧	38.37.04	3,1 MM	13 MM	70 MM
WEDGE Ø4,1	⑩	38.37.05	4,1 MM	13 MM	70 MM
CONIQUE	⑪	38.37.20	3,1 MM	13 MM	70 MM

Available on special order :  
 - Burrs length 50 mm et 65 mm  
 - Tip Stryker

# PERCUTANEOUS INSTRUMENTS



N°	PRODUCT	REF.
1	Elevator	31.33.02
2	Bone lever	31.33.03
3	Bone rasp	31.33.01
4	Surgical beaver blade holder 15.5cm	6061
5	Beaver blades (box of 25) SM64	5904
6	Storage box	14.10.99


# OUR RANGE OF SCREWS

SCREWDRIVER



Ø2.0

Ø2.5



Ø 3.0 Ø 4.0

**MI SCREW**

Ø2.5

Ø4.0



Ø 4.6 Ø 7.0

**EAGLE**

Ø2.5

Ø4.0



Ø 4.6 Ø 7.0

**ARTHROSTAB**

SCREWDRIVER



Ø2.0



Ø 2.0

**MISS**

Ø2.0



Ø2.0

**S.B.S.**

Torx T7



Ø 2.5

**MY SCREW**

Ø2.0



Ø 2.6

**PPDA**

Ø2.0



Ø 3.0

# OUR REFERENCES


OrthoCAPE

SCREWS				
NAME	Ø	LONG.		REF.
MI SCREW Ø3.0	3.0	18 to 60 MM		38.92.xx
MI SCREW Ø4.0	4.0	32 to 60 MM		38.23.xx
EAGLE Ø4.6	4.6	20 to 60 MM		38.98.xx
EAGLE Ø7.0	7.0	30 to 90 MM	THREAD 16	38.99.xxA
EAGLE Ø7.0	7.0	55 to 90 MM	THREAD 22	38.99.xxB
EAGLE Ø7.0	7.0	70 to 90 MM	THREAD 28	38.99.xxC
ARTHROSTAB Ø4.6	4.6	20 to 60 MM		38.21.xx
ARTHROSTAB Ø7.0	7.0	30 to 90 MM		38.22.xx
MISS	2.0	13 to 22 MM		38.52.xx
S.B.S.	2.0	5 to 24 MM		14.52.xx
MY SCREW Ø2.5	2.5	12 to 34 MM		38.88.xx
MY SCREW Ø3.0	3.0	10 to 34 MM		38.27.xx
PPDA Ø2.6	2.6	12 to 34 MM		38.49.xx
PPDA Ø3.0	3.0	12 to 34 MM		38.65.xx
Pin	1.0	70 MM	2 TIPS	14.63.01
Pin	1.0	150 MM	1 TIP	14.63.48
Pin	1.50	150 MM	1 TIP	14.63.44
Pin	2.2	200 MM	1 TIP	14.63.10

# MI SCREW

Ø 3.0 MM Ø 4.0 MM


- Standard ISO 5832-3 Titanium screw
- Self-drilling self-tapping screw
- Forefoot
- Sterile packaging
- Pin (K-Wire) Ø 1.0 for Ø 3.0 mm screw
- Pin (K-Wire) Ø 1.5 for Ø 4.0 mm screw



	Ø 3.0	Ø 4.0
LONG.	REF.	REF.
18	38.92.18	
20	38.92.20	
22	38.92.22	
24	38.92.24	
26	38.92.26	
28	38.92.28	
30	38.92.30	
32	38.92.32	38.23.32
34	38.92.34	38.23.34
36	38.92.36	38.23.36
38	38.92.38	38.23.38
40	38.92.40	38.23.40
42	38.92.42	38.23.42
44	38.92.44	38.23.44
46	38.92.46	38.23.46
48	38.92.48	38.23.48
50	38.92.50	38.23.50
52	38.92.52	38.23.52
54	38.92.54	38.23.54
56	38.92.56	38.23.56
58	38.92.58	38.23.58
60	38.92.60	38.23.60

## INSTRUMENTS


N°	PRODUCTS	REF.
1	Drill head for Ø 4.0 screw	14.37.21
2	Screwdriver blade for Ø 4.0 screw	14.30.10
3	Drill for Ø 4.0 screw	14.37.33
4	Screwdriver handle AO	14.32.03
5	K-wire Ø 1.50 Lg 150 mm 1 tip	14.63.44
6	Screwdriver blade for Ø 3.0 screw	14.30.23
7	Drill for Ø 3.0 screw	14.37.24
8	Measurer	14.33.33
9	Drill for head for Ø 3.0 screw	14.37.20
10	K-wire Ø 1.0 Lg 150 mm 1 tip	14.63.48
11	Storage box	14.10.01



# EAGLE

Ø 4.6 MM Ø 7.0 MM

- Standard ISO 5832-3 Titanium screw
- Mid and hindfoot
- Sterile packaging
- Pin (K-Wire) Ø 1.5 for Ø 4.6 mm screw  
Pin (K-Wire) Ø 2.2 for Ø 7.0 mm screw



Ø 4.6 mm

Ø 7.0 mm

Ø 4.6 mm	
LONG.	REF.
20	38.98.20
22	38.98.22
24	38.98.24
26	38.98.26
28	38.98.28
30	38.98.30
32	38.98.32
34	38.98.34
36	38.98.36
38	38.98.38
40	38.98.40
42	38.98.42
44	38.98.44
46	38.98.46
48	38.98.48
50	38.98.50
52	38.98.52
54	38.98.54
56	38.98.56
58	38.98.58
60	38.98.60

Ø 7.0 mm					
LONG.	REF.	LONG.	REF.	LONG.	REF.
THREAD 16 MM		THREAD 22 MM		THREAD 28 MM	
30	38.99.30A				
35	38.99.35A				
40	38.99.40A				
45	38.99.45A				
50	38.99.50A				
55	38.99.55A	55	38.99.55B		
60	38.99.60A	60	38.99.60B		
65	38.99.65A	65	38.99.65B		
70	38.99.70A	70	38.99.70B	70	38.99.70C
75	38.99.75A	75	38.99.75B	75	38.99.75C
80	38.99.80A	80	38.99.80B	80	38.99.80C
85	38.99.85A	85	38.99.85B	85	38.99.85C
90	38.99.90A	90	38.99.90B	90	38.99.90C

# EAGLE

Ø 4.6 MM Ø 7.0 MM




## Ø 4.6 MM

### INSTRUMENTS

PRODUCTS	REF.
1 Screwdriver handle AO	14.32.03
2 Screwdriver blade for Ø 4.6 screw	14.30.17
3 Drill head for Ø 4.6 screw	14.37.34
4 Tap Eagle Arthrostab Ø 4.6 screw	14.33.36
5 Drill for Ø 4.6 screw	14.37.33
6 Measurer	14.33.33
7 Pins Ø 1.5 Lg 150 mm 1 punta trocar	14.63.44
8 Storage box	14.10.96

screws are used with 1 tip pin  
Ø 1.5 Lg 150 mm



## Ø 7.0 MM

### INSTRUMENTS


PRODUCTS	REF.
1 Screwdriver handle Eagle Arthrostab Ø 7.0 screw	14.32.04
2 Drill head for Ø 7.0 screw	14.37.35
3 Screwdriver blade canulated for Ø 7.0 screw	14.30.18
4 K-wire Ø 2.2 Lg 200 mm 1 tip	14.63.10
5 Measurer	14.33.34
6 Canulated drill for Ø 7.0 screw	14.37.36
7 Tap Eagle Arthrostab Ø 7.0 screw	14.33.35
8 Screw driver blade for Ø 7.0 screw	14.30.19
9 Storage box	14.10.97

screws are used with 1 tip pin  
Ø 2.2 Lg 200 mm

# ARTHROSTAB

Ø 4.6 MM    Ø 7.0 MM

- Standard ISO 5832-3 Titanium screw
- Mid and hindfoot
- Sterile packaging
- Pin (K-Wire) Ø 1.5 for Ø 4.6 mm screw  
Pin (K-Wire) Ø 2.2 for Ø 7.0 mm screw



Ø 4.6 mm

Ø 7.0 mm

## Ø 4.6 mm

LONG.	REF.
20	38.21.20
22	38.21.22
24	38.21.24
26	38.21.26
28	38.21.28
30	38.21.30
32	38.21.32
34	38.21.34
36	38.21.36
38	38.21.38
40	38.21.40
42	38.21.42
44	38.21.44
46	38.21.46
48	38.21.48
50	38.21.50
52	38.21.52
54	38.21.54
56	38.21.56
58	38.21.58
60	38.21.60

## Ø 7.0 mm

LONG.	REF.
30	38.22.30
35	38.22.35
40	38.22.40
45	38.22.45
50	38.22.50
55	38.22.55
60	38.22.60
65	38.22.65
70	38.22.70
75	38.22.75
80	38.22.80
85	38.22.85
90	38.22.90

# ARTHROSTAB

Ø 4.6 MM Ø 7.0 MM




## Ø 4.6 MM

### INSTRUMENTS

	PRODUCTS	REF.
1	Screwdriver handle AO	14.32.03
2	Screwdriver blade for Ø 4.6 screw	14.30.17
3	Drill head for Ø 4.6 screw	14.37.34
4	Tap Eagle Arthrostab Ø 4.6 screw	14.33.36
5	Drill for Ø 4.6 screw	14.37.33
6	Measurer	14.33.33
7	Pins Ø 1.5 Lg 150 mm 1 punta trocar	14.63.44
8	Storage box	14.10.96

screws are used with 1 tip pin

Ø 1.5 Lg 150 mm



## Ø 7.0 MM

### INSTRUMENTS

	PRODUCTS	REF.
1	Screwdriver handle Eagle Arthrostab Ø 7.0 screw	14.32.04
2	Drill head for Ø 7.0 screw	14.37.35
3	Screwdriver blade canulated for Ø 7.0 screw	14.30.18
4	K-wire Ø 2.2 Lg 200 mm 1 tip	14.63.10
5	Measurer	14.33.34
6	Canulated drill for Ø 7.0 screw	14.37.36
7	Tap Eagle Arthrostab Ø 7.0 screw	14.33.35
8	Screw driver blade for Ø 7.0 screw	14.30.19
9	Storage box	14.10.97

screws are used with 1 tip pin


Ø 2.2 Lg 200 mm

# MY SCREW

Ø 2.5 MM

Head Ø 3.5 mm


- Standard ISO 5832-3 Titanium screw
- Self-drilling self-tapping screw
- Forefoot
- Sterile packaging
- K-wire Ø 1.0 for Ø 2.5 & 3.0 mm screw



	Ø2.5
LONG.	REF.
10	38.88.10
12	38.88.12
14	38.88.14
16	38.88.16
18	38.88.18
20	38.88.20
22	38.88.22
24	38.88.24
26	38.88.26
28	38.88.28
30	38.88.30
32	38.88.32
34	38.88.34


## INSTRUMENTS

N°	PRODUCTS	REF.
1	Canulated drill for Ø2.5 screw	14.37.25
2	Canulated drill XS for Ø2.5 screw	14.37.27
3	Measurer	14.33.22
4	Screwdriver blade Torx 7	14.30.22
5	Screwdriver blade Torx 7 canulated	14.30.21
6	K-wire Ø 1.0 Lg 70 mm	14.63.01
7	Screwdriver handle AO	14.32.03
8	Storage box	14.10.88



# PPDA

Ø 2.6 MM    Ø 3.0 MM



Ø 2.6 mm


Ø 3.0 mm

- Standard ISO 5832-3 Titanium screw
- The Self-drilling PPDA® Screw is a Perforated Auto Compression Variable Pitch screw
- Forefoot
- Sterile packaging
- K-wire Ø 1.0 for Ø 2.6 et Ø 3.0 screw

Ø 2.6 mm		Ø 3.0 mm	
LONG.	REF.	LONG.	REF.
12	38.49.12	12	38.65.12
14	38.49.14	14	38.65.14
16	38.49.16	16	38.65.16
18	38.49.18	18	38.65.18
20	38.49.20	20	38.65.20
22	38.49.22	22	38.65.22
24	38.49.24	24	38.65.24
26	38.49.26	26	38.65.26
28	38.49.28	28	38.65.28
30	38.49.30	30	38.65.30
32	38.49.32	32	38.65.32
34	38.49.34	34	38.65.34

## INSTRUMENTS

N°	PRODUCTS	REF.
1	Canulated drill for PPDA Ø 2.6	14.37.25
2	Canulated drill XS for PPDA Ø 2.6	14.37.27
3	Measurer	14.33.22
4	Screwdriver blade for PPDA Ø2.6 & Ø3.0	14.30.23
5	Canulated drill for PPDA Ø 3.0	14.37.23
6	K-wire Ø 1.0 Lg 70 mm	14.63.01
7	Screwdriver handle AO	14.32.03
8	Storage box	14.10.88



Version B

# MISS

Ø 2.0 MM

- Standard ISO 5832-3 Titanium screw
- Twist-off, self-drilling and self-tapping screw
- Forefoot
- Sterile packaging



Ø 2.0 mm

Ø 2.0	
LONG.	REF.
13	38.52.13
15	38.52.15
17	38.52.17
19	38.52.19
22	38.52.22

## INSTRUMENTS

N°	PRODUCTS	REF.
1	Trocar MISS	14.33.25
2	Screwdriver blade MISS	14.30.12
3	Screwdriver handle MISS	14.31.12
4	Cutting guide left	14.33.39
5	Cutting guide right	14.33.39



1

**Trocar MISS**

14.33.25



2

**Screwdriver blade MISS**

14.30.12

3



**Screwdriver handle AO**

4

**Cutting guide left**

14.33.39

5

**Cutting guide right**

14.33.39



# SBS

Ø 2.0

- Standard ISO 5832-3 Titanium screw
- Twist-off, self-drilling and self-tapping screw
- Forefoot
- Sterile packaging



Ø 2.0mm

Ø 2.0 mm	
LONG.	RÉF.
05	14.52.05
06	14.52.06
07	14.52.07
08	14.52.08
09	14.52.09
10	14.52.10
11	14.52.11
12	14.52.12
13	14.52.13
14	14.52.14
15	14.52.15
16	14.52.16
17	14.52.17
18	14.52.18
19	14.52.19
20	14.52.20
21	14.52.21
22	14.52.22
23	14.52.23
24	14.52.24

## INSTRUMENTS

N°	PRODUCTS	REF.
1	Screwdriver blade SBS	14.30.02
2	Screwdriver handle AO	14.32.03



# THREADED PIN

LENGHT 150 MM

- 316L ISO 5832-1 pin
- Extremitie triangular
- Singly sterile packaging Pin + ball



THREADED PINS AND ANTI-MIGRATION BALL

**LENGHT. 150 mm**

DIAMETER	REF.
Threaded pins Ø 1.0 mm + anti-migration ball	38.63.13/18
Threaded pins Ø 1.2 mm + anti-migration ball	38.63.14/18
Threaded pins Ø 1.5 mm + anti-migration ball	38.63.15/19
Threaded pins Ø 1.8 mm + anti-migration ball	38.63.16/20
Threaded pins Ø 2.0 mm + anti-migration ball	38.63.17/20

# VARISATION STAPLES

ANGLE 26° OR 90°

- 316L ISO 5832-1 stainless steel
- Sterile packaging
- Osteotomy AKIN

	ANGLE 26°	ANGLE 90°
DISTANCE 8 mm	14.55.04	14.55.02
DISTANCE 10 mm	14.55.03	14.55.01




ANGLE 26°



ANGLE 90°

## INSTRUMENTS

N°	PRODUCTS	REF.
1	Drilling guide	14.33.55
2	Positioner 26°	14.33.52
3	Impactor 26°	14.33.54
4	Positioner 90°	14.33.51
5	Impactor 90°	14.33.53
6	Storage box	14.10.50




# JACK-S SCREW

6 SIZES

- Titanium screw
- Subtalar screw for calcaneal repositioning
- Easy insertion / Excellent stability
- Treatment of child and adult flatfoot



REF	Ø	LONG.
42-0712	7.0	12
42-0814	8.0	14
42-0914	9.0	14
42-1014	10.0	14
42-1116	11.0	16
42-1216	12.0	16



Download CE certificate



**[www.ortho-cape.com](http://www.ortho-cape.com)**



IBAN FR76 1027 8022 0000 0210 4320 276 - BIC CMCIFR2A  
N°VAT - FR 61 890810542  
SIRET 890 810 542 00021 - APE