


MY SCREW PLEINE

Ø 2.5 MM NON CANULÉ

Tête Ø 3.5 mm



Ø 2.5 mm

- Vis en titane ISO 5832-3
- Auto taraudeuse, auto foreuse et auto compressive
- Avant-pied
- Conditionnement stérile

Ø 2.5 mm	
LONG.	RÉF.
10	38.24.10
12	38.24.12
14	38.24.14
16	38.24.16
18	38.24.18
20	38.24.20
22	38.24.22
24	38.24.24
26	38.24.26
28	38.24.28
30	38.24.30
32	38.24.32
34	38.24.34

ANCILLAIRE

N°	PRODUIT	RÉF.
1	Forêt pour vis My Screw Ø2.5	14.37.25
2	Forêt XS pour vis My Screw Ø2.5	14.37.27
3	Tournevis (lame) 6 pans	14.30.23A
4	Manche de tournevis AO universel	14.32.03
5	Boite de rangement	14.10.88

