


# MY SCREW

Ø 2.5 MM Ø3.0

Tête Ø 3.5 mm

Tête Ø 4.0 mm



Ø 2.5 mm

Ø 3.0 mm

- Vis en titane ISO 5832-3
- Auto taraudeuse, auto foreuse et auto compressive
- Avant-pied
- Conditionnement stérile
- Broche Ø 1.0 pour vis Ø 2.5 mm et 3.0 mm

Ø 2.5 mm		Ø 3.0 mm	
LONG.	RÉF.	LONG.	RÉF.
12	38.88.12		
14	38.88.14	14	38.27.14
16	38.88.16	16	38.27.16
18	38.88.18	18	38.27.18
20	38.88.20	20	38.27.20
22	38.88.22	22	38.27.22
24	38.88.24	24	38.27.24
26	38.88.26	26	38.27.26
28	38.88.28	28	38.27.28
30	38.88.30	30	38.27.30
32	38.88.32	32	38.27.32
34	38.88.34	34	38.27.34

## ANCILLAIRE

N°	PRODUIT	RÉF.
1	Forêt canulé pour vis My Screw	14.37.25
2	Forêt canulé XS pour vis My Screw	14.37.27
3	Mesureur sur broche 70 mm	14.33.22
4	Tournevis (lame) Torx 7 canulé pour Ø2.5	14.30.21
5	Tournevis (lame) Torx 8 canulé pour Ø3.0	14.30.27
6	Broches Ø 1.0 Lg 70 mm	14.63.01
7	Manche de tournevis AO universel	14.32.03
8	Boite de rangement	14.10.88

