

SBS

Ø 2.0

■ Vis en titane ISO 5832-3

■ Sécable, auto perforante et auto taraudeuse

■ Avant-pied

■ Conditionnement stérile

Tête Ø 3.0 mm



Ø 2.0mm



Ø 2.0 mm	
LONG.	RÉF.
09	14.52.09
10	14.52.10
11	14.52.11
12	14.52.12
13	14.52.13
14	14.52.14
15	14.52.15
16	14.52.16
17	14.52.17
18	14.52.18
19	14.52.19

ANCILLAIRE

N°	PRODUIT	RÉF.
1	Manche de tournevis AO	14.32.03
2	Lame de tournevis S.B.S	14.30.02
3	Boite	14.10.90

