



Ortho CAPE


FOOT SURGERY

For any questions regarding our products,
contact us at info@ortho-cape.com

PERCUTANEOUS BURRS

OrthoCAPE

- X90CrMoV18 Steel burrs
- Percutaneous surgery
- Singly sterile packaging
- Standards burrs Ø 2.32 Lg 70



	N°	REF.	Ø	CUTTING LG	TOTAL LG
SHANNON SHORT	①	38.37.01	2 MM	8 MM	70 MM
SHANNON LENG	②	38.37.02	2 MM	12 MM	70 MM
SHANNON ISHAM	③	38.37.03	2 MM	12 MM	70 MM
SHANNON XL	④	38.37.07	2 MM	22 MM	70 MM
SHANNON XXL	⑤	38.37.08	3 MM	20 MM	70 MM
SHANNON XXL LG 30	⑥	38.37.09	3 MM	30 MM	100 MM
BROPHY	⑦	38.37.06	5 MM	15 MM	70 MM
BROPHY SHORT	⑧	38.37.0X	5 MM	13 MM	70 MM
BROPHY SLIM	⑨	38.37.0X	4 MM	13 MM	70 MM
WEDGE Ø3,1	⑩	38.37.04	3,1 MM	13 MM	70 MM
WEDGE Ø4,1	⑪	38.37.05	4,1 MM	13 MM	70 MM
CONIQUE	⑫	38.37.20	3,1 MM	13 MM	70 MM
BALL Ø4,0	⑬	38.37.XX	4 MM	3 MM	70 MM
BALL Ø5,0	⑭	38.37.XX	5 MM	4 MM	70 MM

coming soon
coming soon

coming soon
coming soon

PERCUTANEOUS INSTRUMENTS

OrthoCAPE



N°	PRODUCT	REF.
1	Elevator	31.33.02
2	Bone lever	31.33.03
3	Bone rasp	31.33.01
4	Surgical beaver blade holder 15.5cm	6061
5	Beaver blades (box of 25) SM64	5904
6	Storage box	14.10.99

OUR RANGE OF SCREWS

OrthoCAPE

SCREWDRIVER

Ø2.0

Ø2.5

Ø2.5

Ø4.0


Ø1.6

T7

T8


Ø2.5

Ø4.0




Ø 3.0 Ø 4.0

MI SCREW




Ø 4.6 Ø 7.0

EAGLE



Ø 2.5 Ø 2.5 Ø 3.0

MINI ARTHROSTAB



Ø 4.6 Ø 7.0

ARTHROSTAB

SCREWDRIVER

Ø2.0

Ø2.0

Ø1.6

T7

T8

Ø2.0

Ø2.0




Ø 2.0

MISS




Ø 2.0

SBS



Ø 2.5 Ø 2.5 Ø 3.0

MY SCREW



Ø 2.6 Ø 3.0

PPDA

OUR SCREWS REFERENCES


OrthoCAPE

SCREWS				
NAME	Ø	LENG.		REF.
MI SCREW Ø3.0	3.0	18 to 60 MM		38.92.xx
MI SCREW Ø4.0	4.0	32 to 60 MM		38.23.xx
EAGLE Ø4.6	4.6	20 to 60 MM		38.98.xx
EAGLE Ø7.0	7.0	30 to 90 MM	THREAD 16	38.99.xxA
EAGLE Ø7.0	7.0	55 to 90 MM	THREAD 22	38.99.xxB
EAGLE Ø7.0	7.0	70 to 90 MM	THREAD 28	38.99.xxC
MINI ARTHROSTAB Ø2.5	2.5	12 to 34 MM		38.25.xx
MINI ARTHROSTAB Ø3.0	3.0	14 to 40 MM		38.26.xx
BEVELED MINI ARTHROSTAB	2.5	14 to 34 MM		38.28.xx
ARTHROSTAB Ø4.6	4.6	20 to 60 MM		38.21.xx
ARTHROSTAB Ø7.0	7.0	30 to 90 MM		38.22.xx
MISS	2.0	13 to 22 MM		38.52.xx
S.B.S.	2.0	5 to 24 MM		14.52.xx
MY SCREW	2.5	12 to 34 MM		38.88.xx
MY SCREW US	2.5	10 to 34 MM		38.24.xx
MY SCREW	3.0	14 to 34 MM		38.27.xx
PPDA Ø2.6	2.6	12 to 34 MM		38.49.xx
PPDA Ø3.0	3.0	12 to 34 MM		38.65.xx
Pin	0.8	70 MM		14.63.41
Pin	1.0	70 MM		14.63.01
Pin	1.0	150 MM	1 TIP	14.63.48
Pin	1.50	150 MM	1 TIP	14.63.44
Pin	2.2	200 MM	1 TIP	14.63.10

MI SCREW

Ø 3.0 MM Ø 4.0 MM

- Standard ISO 5832-3 Titanium screw
- Self-drilling self-tapping screw
- Forefoot
- Sterile packaging
- Ø 1.0 Pin for Ø 3.0 mm screw
- Ø 1.5 Pin for Ø 4.0 mm screw




Ø 3.0 mm Ø 4.0 mm

	Ø 3.0	Ø 4.0
LENG.	REF.	REF.
18	38.92.18	
20	38.92.20	
22	38.92.22	
24	38.92.24	
26	38.92.26	
28	38.92.28	
30	38.92.30	
32	38.92.32	38.23.32
34	38.92.34	38.23.34
36	38.92.36	38.23.36
38	38.92.38	38.23.38
40	38.92.40	38.23.40
42	38.92.42	38.23.42
44	38.92.44	38.23.44
46	38.92.46	38.23.46
48	38.92.48	38.23.48
50	38.92.50	38.23.50
52	38.92.52	38.23.52
54	38.92.54	38.23.54
56	38.92.56	38.23.56
58	38.92.58	38.23.58
60	38.92.60	38.23.60

INSTRUMENTS

N°	PRODUCTS	REF.
1	Ø 4.0 screw Drill head	14.37.21
2	Ø 4.0 screw Screwdriver blade	14.30.10
3	Ø 4.0 screw Drill	14.37.33
4	AO Screwdriver handle	14.32.03
5	Ø 1.50 Lg 150 mm 1 tip Pin	14.63.44
6	Ø 3.0 screw Screwdriver blade	14.30.23
7	Ø 3.0 screw Drill	14.37.24
8	Measurer	14.33.33
9	Ø 3.0 screw Drill for head	14.37.20
10	Ø 1.0 Lg 150 mm 1 tip Pin	14.63.48
11	Storage box	14.10.01



EAGLE

Ø 4.6 MM Ø 7.0 MM

- Standard ISO 5832-3 Titanium screw
- Mid and hindfoot
- Sterile packaging
- Ø 1.5 Pin for Ø 4.6 mm screw
- Ø 2.2 Pin for Ø 7.0 mm screw

Head Ø 6.0 mm



Head Ø 9.5 mm



Ø 4.6 mm	
LENG.	REF.
20	38.98.20
22	38.98.22
24	38.98.24
26	38.98.26
28	38.98.28
30	38.98.30
32	38.98.32
34	38.98.34
36	38.98.36
38	38.98.38
40	38.98.40
42	38.98.42
44	38.98.44
46	38.98.46
48	38.98.48
50	38.98.50
52	38.98.52
54	38.98.54
56	38.98.56
58	38.98.58
60	38.98.60

Ø 7.0 mm					
LENG.	REF.	LENG.	REF.	LENG.	REF.
THREAD 16 MM		THREAD 22 MM		THREAD 28 MM	
30	38.99.30A				
35	38.99.35A				
40	38.99.40A				
45	38.99.45A				
50	38.99.50A				
55	38.99.55A	55	38.99.55B		
60	38.99.60A	60	38.99.60B		
65	38.99.65A	65	38.99.65B		
70	38.99.70A	70	38.99.70B	70	38.99.70C
75	38.99.75A	75	38.99.75B	75	38.99.75C
80	38.99.80A	80	38.99.80B	80	38.99.80C
85	38.99.85A	85	38.99.85B	85	38.99.85C
90	38.99.90A	90	38.99.90B	90	38.99.90C

Ø 4.6 mm

Ø 7.0 mm

EAGLE

Ø 4.6 MM Ø 7.0 MM



Ø 4.6 MM

INSTRUMENTS

	PRODUCTS	REF.
1	AO Screwdriver handle	14.32.03
2	Ø 4.6 screw Screwdriver blade	14.30.17
3	Ø 4.6 screw Drill head	14.37.34
4	Ø 4.6 screw Tap	14.33.36
5	Ø 4.6 screw Drill	14.37.33
6	Measurer	14.33.33
7	Ø 1.5 Lg 150 mm 1 trocar tip Pin	14.63.44
8	Storage box	14.10.96

screws are used with 1 tip pin
Ø 1.5 Lg 150 mm



Ø 7.0 MM

INSTRUMENTS


	PRODUCTS	REF.
1	Ø 7.0 screw Screwdriver handle	14.32.04
2	Ø 7.0 screw Drill head for	14.37.35
3	Canulated Ø 7.0 screw Screwdriver blade	14.30.18
4	Ø 2.2 Lg 200 mm 1 tip Pin	14.63.10
5	Measurer	14.33.34
6	Ø 7.0 screw Canulated drill	14.37.36
7	Ø 7.0 screw Tap	14.33.35
8	Ø 7.0 screw Screwdriver blade	14.30.19
9	Storage box	14.10.97

screws are used with 1 tip pin
Ø 2.2 Lg 200 mm

ARTHROSTAB

Ø 4.6 MM Ø 7.0 MM

- Standard ISO 5832-3 Titanium screw
- Mid and hindfoot
- Sterile packaging
- Ø 1.5 Pin for Ø 4.6 mm screw
- Ø 2.2 Pin for Ø 7.0 mm screw



Ø 4.6 mm	
LENG.	REF.
20	38.21.20
22	38.21.22
24	38.21.24
26	38.21.26
28	38.21.28
30	38.21.30
32	38.21.32
34	38.21.34
36	38.21.36
38	38.21.38
40	38.21.40
42	38.21.42
44	38.21.44
46	38.21.46
48	38.21.48
50	38.21.50
52	38.21.52
54	38.21.54
56	38.21.56
58	38.21.58
60	38.21.60

Ø 7.0 mm	
LENG.	REF.
30	38.22.30
35	38.22.35
40	38.22.40
45	38.22.45
50	38.22.50
55	38.22.55
60	38.22.60
65	38.22.65
70	38.22.70
75	38.22.75
80	38.22.80
85	38.22.85
90	38.22.90

ARTHROSTAB

Ø 4.6 MM Ø 7.0 MM



Ø 4.6 MM

INSTRUMENTS

	PRODUCTS	REF.
1	AO Screwdriver handle	14.32.03
2	Ø 4.6 screw Screwdriver blade	14.30.17
3	Ø 4.6 screw Drill head	14.37.34
4	Ø 4.6 screw Tap	14.33.36
5	Ø 4.6 screw Drill	14.37.33
6	Measurer	14.33.33
7	Ø 1.5 Lg 150 mm 1 trocar tip Pin	14.63.44
8	Storage box	14.10.96

screws are used with 1 tip pin
Ø 1.5 Lg 150 mm



Ø 7.0 MM

INSTRUMENTS

	PRODUCTS	REF.
1	Ø 7.0 screw Screwdriver handle	14.32.04
2	Ø 7.0 screw Drill head for	14.37.35
3	Canulated Ø 7.0 screw Screwdriver blade	14.30.18
4	Ø 2.2 Lg 200 mm 1 tip Pin	14.63.10
5	Measurer	14.33.34
6	Ø 7.0 screw Canulated drill	14.37.36
7	Ø 7.0 screw Tap	14.33.35
8	Ø 7.0 screw Screwdriver blade	14.30.19
9	Storage box	14.10.97

screws are used with 1 tip pin
Ø 2.2 Lg 200 mm

MINI ARTHROSTAB

Ø 2.5 MM Ø 3.0 MM

Head Ø 3.5 mm

Head Ø 4.0 mm

Ø 2.5 mm

Ø 3.0 mm


- Standard ISO 5832-3 Titanium screw
- Forefoot
- Sterile packaging
- Ø 1.0 Pin for Ø 2.5 & 3.0 mm screw

Ø 2.5 mm	
LENG.	REF.
12	38.25.12
14	38.25.14
16	38.25.16
18	38.25.18
20	38.25.20
22	38.25.22
24	38.25.24
26	38.25.26
28	38.25.28
30	38.25.30
32	38.25.32
34	38.25.34

Ø 3.0 mm	
LENG.	REF.
14	38.26.14
16	38.26.16
18	38.26.18
20	38.26.20
22	38.26.22
24	38.26.24
26	38.26.26
28	38.26.28
30	38.26.30
32	38.26.32
34	38.26.34


INSTRUMENTS

N°	PRODUCT	REF.
1	Mini arthrostab canulated drill	14.37.25
2	Mini Arthrostab XS canulated drill	14.37.27
3	Measurer 70 mm	14.33.22
4	Ø2.5 canulated Torx 7 screwdriver blade	14.30.21
5	Ø3.0 Torx 8 canulated screwdriver blade	14.30.27
6	Ø 1.0 Lg 70 mm Pin	14.63.01
7	AO Screwdriver handle	14.32.03
8	Storage box	14.10.88



MINI ARTHROSTAB

Ø 2.5 MM BEVELED




- Standard ISO 5832-3 Titanium screw
- Forefoot
- Sterile packaging
- Ø 1.0 Pin for Ø 2.5 mm screw

Ø 2.5 mm	
LENG.	REF.
14	38.28.14
16	38.28.16
18	38.28.18
20	38.28.20
22	38.28.22
24	38.28.24
26	38.28.26
28	38.28.28
30	38.28.30
32	38.28.32
34	38.28.34

INSTRUMENTS

N°	PRODUCT	REF.
1	Mini Arthrostab drill	14.37.25
2	Mini Arthrostab XS drill	14.37.27
3	Mini Arthrostab hex screwdriver blade	14.30.23
4	AO Screwdriver handle	14.32.03
5	Storage box	14.10.88



MISS

Ø 2.0 MM

Head Ø 2.7 mm

- Standard ISO 5832-3 Titanium screw
- Twist-off, self-drilling and self-tapping screw
- Forefoot
- Sterile packaging

Ø 2.0	
LENG.	REF.
13	38.52.13
15	38.52.15
17	38.52.17
19	38.52.19
22	38.52.22

Ø 2.0 mm

INSTRUMENTS

N°	PRODUCT	REF.
1	AO Screwdriver handle	14.32.03
2	MISS AO Screwdriver blade	14.30.12
3	Removal screwdriver blade	14.30.03
4	MISS square-tip screwdriver blade	14.30.25
5	Storage box	14.10.90



SBS

Ø 2.0

- Standard ISO 5832-3 Titanium screw
- Twist-off, self-drilling and self-tapping screw
- Forefoot
- Sterile packaging

Head Ø 3.0 mm



Ø 2.0 mm	
LENG.	REF.
05	14.52.05
06	14.52.06
07	14.52.07
08	14.52.08
09	14.52.09
10	14.52.10
11	14.52.11
12	14.52.12
13	14.52.13
14	14.52.14
15	14.52.15
16	14.52.16
17	14.52.17
18	14.52.18
19	14.52.19
20	14.52.20
21	14.52.21
22	14.52.22
23	14.52.23
24	14.52.24

INSTRUMENTS

N°	PRODUCT	REF.
1	AO Screwdriver handle	14.32.03
2	S.B.S screwdriver blade	14.30.02
3	Storage box	14.10.90




MY SCREW

Ø 2.5 MM

- Standard ISO 5832-3 Titanium screw
- Self-drilling self-tapping screw
- Forefoot
- Sterile packaging
- Ø 1.0 Pin for Ø 2.5 & 3.0 mm screw

Head Ø 4.0 mm

Head Ø 3.5 mm




Ø 2.5 mm	
LENG.	REF.
12	38.88.12
14	38.88.14
16	38.88.16
18	38.88.18
20	38.88.20
22	38.88.22
24	38.88.24
26	38.88.26
28	38.88.28
30	38.88.30
32	38.88.32
34	38.88.34

Ø 3.0 mm	
LENG.	REF.
14	38.27.14
16	38.27.16
18	38.27.18
20	38.27.20
22	38.27.22
24	38.27.24
26	38.27.26
28	38.27.28
30	38.27.30
32	38.27.32
34	38.27.34

INSTRUMENTS


N°	PRODUCTS	REF.
1	My Screw Canulated drill	14.37.25
2	XS My Screw Canulated drill	14.37.27
3	Measurer	14.33.22
4	Ø2.5 Canulated Torx 7 screwdriver blade	14.30.21
5	Ø3.0 Canulated Torx 8 screwdriver blade	14.30.27
6	Ø 1.0 Lg 70 mm Pin	14.63.01
7	AO Screwdriver handle	14.32.03
8	Storage box	14.10.88



MY SCREW US

Ø 2.5 MM UNDRILLED SCREW

Head Ø 3.5 mm




Ø 2.5 mm

- Standard ISO 5832-3 Titanium screw
- Self-drilling self-tapping screw
- Forefoot
- Sterile packaging

Ø 2.5 mm	
LENG.	REF.
10	38.24.10
12	38.24.12
14	38.24.14
16	38.24.16
18	38.24.18
20	38.24.20
22	38.24.22
24	38.24.24
26	38.24.26
28	38.24.28
30	38.24.30
32	38.24.32
34	38.24.34

INSTRUMENTS

N°	PRODUCT	REF.
1	My Screw drill	14.37.25
2	XS My Screw drill	14.37.27
3	My Screw hex screwdriver blade	14.30.23
4	AO Screwdriver handle	14.32.03
5	Storage box	14.10.90




PPDA

Ø 2.6 MM Ø 3.0 MM

Head Ø 3.6 mm

Head Ø 4.0 mm




Ø 2.6 mm Ø 3.0 mm

- Standard ISO 5832-3 Titanium screw
- The Self-drilling PPDA® Screw is a Perforated Auto Compression Variable Pitch screw
- Forefoot
- Sterile packaging
- Ø 1.0 Pin for Ø 2.6 et Ø 3.0 screw

Ø 2.6 mm		Ø 3.0 mm	
LENG.	REF.	LENG.	REF.
12	38.49.12	12	38.65.12
14	38.49.14	14	38.65.14
16	38.49.16	16	38.65.16
18	38.49.18	18	38.65.18
20	38.49.20	20	38.65.20
22	38.49.22	22	38.65.22
24	38.49.24	24	38.65.24
26	38.49.26	26	38.65.26
28	38.49.28	28	38.65.28
30	38.49.30	30	38.65.30
32	38.49.32	32	38.65.32
34	38.49.34	34	38.65.34

INSTRUMENTS


N°	PRODUCTS	REF.
1	PPDA Ø 2.6 Canulated drill	14.37.25
2	XS PPDA Ø 2.6 Canulated drill	14.37.27
3	Measurer	14.33.22
4	Ø2.6 & Ø3.0 Screwdriver blade	14.30.23
5	PPDA Ø 3.0 Canulated drill	14.37.23
6	Ø 1.0 Lg 70 mm Pin	14.63.01
7	AO Screwdriver handle	14.32.03
8	Storage box	14.10.88



VARISATION STAPLES

ANGLE 26° OR 90°

- 316L ISO 5832-1 stainless steel
- Sterile packaging
- Osteotomy AKIN



ANGLE 26°




ANGLE 90°

	ANGLE 26°	ANGLE 90°
DISTANCE 8 mm	14.55.04	14.55.02
DISTANCE 10 mm	14.55.03	14.55.01

INSTRUMENTS

N°	PRODUCTS	REF.
1	Drilling guide	14.33.55
2	Positioner 26°	14.33.52
3	Impactor 26°	14.33.54
4	Positioner 90°	14.33.51
5	Impactor 90°	14.33.53
6	Storage box	14.10.50




THREADED PIN

LENGHT 150 MM

- 316L ISO 5832-1 pin
- Extremitie triangular
- Singly sterile packaging Pin + ball



THREADED PINS AND ANTI-MIGRATION BALL



LENGHT. 150 mm	
DIAMETER	REF.
Threaded pins Ø 1.0 mm + anti-migration ball	38.63.13/18
Threaded pins Ø 1.2 mm + anti-migration ball	38.63.14/18
Threaded pins Ø 1.5 mm + anti-migration ball	38.63.15/19
Threaded pins Ø 1.8 mm + anti-migration ball	38.63.16/20
Threaded pins Ø 2.0 mm + anti-migration ball	38.63.17/20

JACK-S SCREW


6 SIZES



- Titanium screw
- Subtalar screw for calcaneal repositioning
- Easy insertion / Excellent stability
- Treatment of adult and child flatfoot

REF	Ø	LENG.
42-0712	7.0	12
42-0814	8.0	14
42-0914	9.0	14
42-1014	10.0	14
42-1116	11.0	16
42-1216	12.0	16

REF	PRODUCT
ML-2000	Box
ML-1007	Threaded gauge Ø 7.0
ML-1008	Threaded gauge Ø 8.0
ML-1009	Threaded gauge Ø 9.0
ML-1010	Threaded gauge Ø 10.0
ML-1011	Threaded gauge Ø 11.0
ML-1012	Threaded gauge Ø 12.0
42-4230	Threaded gauge
42-4220	Extractor
42-4240	Probe
14.63.18	Pin Ø 1.8 Lg2 00
42-4243	Dilatator
55-2951	Screwdriver handle




Surgical technic




HAMMERFIX

3 TAILLES



- Zeniva® PEEK
- The implant's lug segment is designed to expand into the proximal phalanx.
- Opposing threads allow for additional compression across the joint.
- Accommodates variations in phalanx size.




	S	M	L
SIZE	2.8	3.4	4.0
OVERALL LENGTH	13.5	15.0	15.0
MAJOR Ø SCREW	2.8	3.4	4.0
MAJOR Ø BARB	4.0	4.5	5.5
BARB LENGTH PROXIMAL	5.5	6.8	6.8
THREADED LENGTH MIDDLE	6.5	7.0	7.0
REF	132-28013-S	132-34015-S	132-40015-S

Surgical technic



[More information](#)



INSTRUMENTS	
REF	PRODUCT
132-00100	Storage box
	Cleaning rod
132-00211	Ø1.1 Guidewire holder S
132-00214	Ø1.4 Guidewire holder M
132-00216	Ø1.6 Guidewire holder L
132-00004	Screwdriver handle
132-00005	Compression forceps
132-00028	Drill/Tap S
132-00034	Drill/Tap M
132-00040	Drill/Tap L
132-00001	Screwdriver S
132-00002	Screwdriver M
132-00003	TScrewdriver L

ENDO-CUT

4 MODELS

- Endoscopic decompression of the interdigital nerve
- Morton's neuroma



45° Angled Micro Hook Blade

14.14.03

Hook Blade

14.14.04

Push/Pull Blade

14.14.02

Pulling Blade

14.14.01

Handle

14.32.10

OrthoCAPE

14.32.10

CARTILAGINATOR

Flat
Single/Double Sided



Flat
Saw Tip



Large

Curved (3D)
Concave/Convex



Small

*Procedure Type	Flat Single-Sided	Flat Double-Sided	Large Curved Concave (3D)	Large Curved Convex (3D)	Small Curved Concave (3D)	Small Curved Convex (3D)	Flat Saw Tip
	168-05010	168-05020	170-01250	170-01251	170-02140	170-02141	170-00011
Ankle Fusion	✓	✓	✓	✓			✓
Charcot	✓	✓					✓
Double Arthrodesis	✓		✓	✓	✓	✓	
Lapidus	✓	✓					✓
Midfoot Fusion	✓	✓					✓
MTP Fusion	✓				✓	✓	
NC Fusion	✓	✓	✓	✓			
Subtalar Fusion	✓	✓	✓	✓			✓
Triple Arthrodesis	✓				✓	✓	✓


More information :
[ortho cape cartilaginator](#)



TENDON REPAIR

HS FIBER®

HS 530



- USP 5 / SIZE 7
- Ø 0.7
- Length : 90 cm
- Color : white & green
- Needles 1/2c 26 mm and 1/2c 48 mm
- Non-absorbable surgical suture UHMWPE
- Double needle sutures

KNEE / ANKLE / HIP



Very easy clamping

USES

- reinforcement and reinsertion ligament and tendinous
- ideal suture and anchorage on knees hip and ankle

Download CE certificate



NEO HS FIBER® LOOP

HS142




- Thread USP 2 : 0,50-0,59 mm
- Length : 90 cm
- Color : white & blue
- Non-absorbable UHMWPE sutures
- Sutures with one 60 mm straight needle
- Loop type: SPLIT-STITCH™

KNEE / HIP



USES

- Grafting, ACL repair, rotator cuff repair, and other soft tissue applications



NEO LOOP 2.0

HS343

- Thread Ø 0,30 mm (USP2/0)
- Length : 60 cm
- Loop length : 5 cm
- Color : black
- Non-absorbable UHMWPE sutures
- Aiguille droite alliage Acier Inox siliconée
- Loop type : SPLIT-STICH™




ANKLE



USES

- Grafting, Achilles tendon repair, and other soft tissue applications



NEO LOOP TAPE 1.4

HS182




- Tape : largeur 1,4 mm
- Length : 60 cm
- Loop length : 5 cm
- Color : black
- Non-absorbable UHMWPE sutures
- Black straight needle made of siliconized stainless-steel alloy
- Loop type : SPLIT-STICH™

KNEE / ANKLE



USES

- Grafting, ACL repair, Achilles tendon repair, and other soft tissue applications



NEO LOOP TAPE 0

HS345



- Thread USP 0
- Length : 60 cm
- Loop length : 5 cm
- Color : black
- Non-absorbable UHMWPE sutures
- Black straight needle made of siliconized stainless-steel alloy
- Loop type : SPLIT-STICH™



KNEE




USES

- Grafting, ACL repair and other soft tissue applications

SHOULDER REPAIR FIBER TAPE HS SUTURE TAPE®

HS 1071



- USP 2 (SIZE 5) 0.50mm to 0.599mm
- 0.50 mm flat ribbon
- Length : 100 cm
- Color: blue
- Needles : 1/2 26 mm and 3/8 40 mm
- Non-absorbable surgical suture tape UHMWPE
- Double needle sutures

SHOULDER



Reinforcement suture in flat ribbon form

Provides a support surface more important while limiting shear risks.

USES

- Ligament reinforcement
- Tendon anchorage
- Reconstruction
- Tissue approximation

Download CE certificate





THE REFERENCE THREAD
FOR
SHOULDER SURGERY



SHOULDER AND ANKLE REPAIR FIBER TAPE

HS SUTURE TAPE®

HS 1072



- USP 2 (SIZE 5) 0.50mm to 0.599mm
- 0.50 mm flat ribbon
- Length : 100 cm
- Color: blue & white
- Needles : 1/2 26 mm and 3/8 40 mm
- Non-absorbable surgical suture tape UHMWPE
- Double needle sutures

SHOULDER / ANKLE



Reinforcement suture in flat ribbon form

Provides a support surface more important while limiting shear risks

Two-tone thread for better visibility under arthroscopy

USES

- Ligament reinforcement
- Tendon anchorage
- Reconstruction
- Tissue approximation

Download CE certificate





TEL. +33 5 63 27 02 25



info@ortho-cape.com



172 rue Lépine
82000 Montauban
FRANCE

Download CE certificate



www.ortho-cape.com

IBAN FR76 1027 8022 0000 0210 4320 276 - BIC CMCIFR2A

N°VAT - FR 61 890810542

SIRET 890 810 542 00021 - APE